

## Blue-Native PAGE

(ref Shagger, H., and Von Jagow, G., Analytical Biochemistry, 1991, 199, 223-231)

### Solutions:

**3xGB** - 3xGel buffer (1.5M aminocaproic acid, 150 mM Bis-tris, pH 7.0)

*19.68g Aminocaproic acid*

*3.14g Bis-tris*

*ad 100 mL H<sub>2</sub>O*

**Cathode buffer** (15mM Bis-tris, 50 mM Tricine, pH 7.0)

*3.14g Bis-tris*

*8.96g Tricine*

*ad 1000mL H<sub>2</sub>O*

**Blue cathode buffer** (Cathode buffer cont. 0.02% Serva Blue G (SBG))

*100 mL Cathode buffer*

*0.02g SBG*

**Anode buffer** (50mM Bis-tris, pH7.0)

*20.93 g Bis-tris*

*ad 2000 mL H<sub>2</sub>O*

**AB** - Acrylamide/bisacrylamide mix (48% acrylamide, 1.5% bisacrylamide (99.5T, C))

*24.0 g acrylamide*

*0.75 g bisacrylamide*

*ad 50 mL H<sub>2</sub>O*

**MB2** - 0.5 ml 3xGB, 0.5 ml 2M Aminocaproic acid, 4µl 500mM EDTA

**SBG** - sample buffer (750 mM Aminocaproic acid, 5% SBG)

*3.75 mL 2M Aminocaproic acid*

*0.5g SBG*

*ad 10 mL H<sub>2</sub>O*

**AC** - 2M aminocaproic acid (13.12g/50 mL)

**LM** - 10% Lauroyl maltoside (0.1g/mL)

**AP** - 10% Amonium persulfate (0.1g/mL)

	Stacking 4% gel		Separating 6% gel		Separating 15% gel	
	2.5 ml	5 ml	5 ml	10 ml	5 ml	10 ml
<b>3xGB</b>	0.82 ml	1.64ml	1.65 ml	3.3 ml	1.65 ml	3.3 ml
<b>AB</b>	0.2 ml	0.4ml	0.6 ml	1.2 ml	1.5 ml	3.0 ml
<b>water</b>	1.45 ml	2.87 ml	2.72 ml	5.44 ml	0.84 ml	1.68 ml
<b>glycerol</b>	0	0	0	0	1 ml	2 ml
<b>AP</b>	30 µl	60 µl	30 µl	60 µl	5 µl	10 µl
<b>TEMED</b>	3 µl	6 µl	2 µl	4 µl	2 µl	2 µl

### Preparation of samples

Use isolated mitochondria or 15 000xg sediments of postnuclear supernatants, or digitonin treated fibroblasts. Perform all steps at 0-4°C:

- resuspend mitochondria in MB2 at 1-3 mg protein/ml
- add LM to final concentration 1%, mix well and leave 15 min
- centrifuge 20 min at 20 000xg
- add SBG to supernatant to final concentration of 0.25%
- apply samples to gel

### Preparation of gel and electrophoresis

Use Miniprotean Bio Rad apparatus with 1 mm spacers and combs:

- mount together the glass "sandwich" cassettes
- prepare 5 ml of 6% and 15% separating gel mixture according to table
- fill the tubing and about 1 cm of cassette height with water
- fill the front reservoir of gradient mixer with 2.8 mL 6% mixture and the rear reservoir with 2.3 mL of 15% gel mixture, put the conical insert into rear reservoir
- fill the gel cassette by means of underlaying the gel under water using peristaltic pump at speed 9
- after polymerization wash gel surface with 1x GB
- prepare 2.5 ml of stacking gel mixture according to table, insert comb and fill up to the top with stacking gel mixture
- after polymerization remove comb and wash wells with 1x GB
- fill the spaces in wells with blue cathode buffer and underlay the samples into wells (sample volume should not exceed 20 µl/well (10 well comb), the amount of protein should be between 10-30 µg prot.)
- place gel cassettes into apparatus and fill upper and lower reservoir with blue cathode and anode buffer
- run 15 min at 40V, then at 80V until the dye reaches 2/3 of the gel, than replace blue cathode with cathode buffer and continue electrophoresis until the dye front reaches the end
- remove gel from between the plates
- use for Western blotting or stain proteins in gel (30 min) and destain the background (1-2h, 2x change of the destain solution)
- for second dimension electrophoresis, cut the gel slice with separated proteins and use directly or freeze and store at - 20°C