

Preparation of cells for Blue Native electrophoresis

- (1) trypsinize cells, wash once with PBS
- (2) resuspend cells in cold PBS (0.5 ml PBS per 1 plate 10cm, 1 ml PBS - 15 cm)
- (3) measure the protein concentration (use Bradford method, sonicate cells)
- (4) add more PBS, spin the cells
- (5) resuspend cells in PBS to final concentration of 5 mg/ml (according the protein concentration measured above) - put them in eppendorf tube
- (6) add equal volume of digitonin (4mg/ml) and incubate on ice for 5 min (final concentration of cells is 2.5mg/ml, final concentration of digitonin is 2 mg/ml, digitonin/protein ratio is 0.8)
- (7) add PBS to final volume app. 1.5 ml
- (8) centrifuge 10 000g, 10 min
- (9) resuspend sediment (membranes) in MB2 buffer (half volume of PBS in step 5)
- (10) add lauryl malthoside (LM, dodecyl malthoside 10%) to final concentration of 1 % (1/10 volume of MB2), incubate on ice 15 min
- (11) centrifuge 20 000g 20min
- (12) put supernatant in a new tube and measure the protein concentration
- (13) add SBG - 1/2 volume of lauryl malthoside
- (14) freeze -20C and run the gel