

# COX ASSAY

## Extraction

Wash cells with cold PBS x2. Place 1.0mL cold PBS on plate and scrape cells. Place in eppendorf tube and centrifuge at 10,000g at 4°C for one minute. Remove supernatant and resuspend pellet in extraction buffer. For a confluent 100mm plate use ~100-200µL. Keep on ice. Homogenize cells with motorized pestle for 15 seconds. Use homogenate within the next few hours. Citrate synthase appears more sensitive to time so perform this assay first.

## Extraction buffer TRA-EDTA pH7.4:

Triethanolamine 50mM (9.28g/L)

EDTA 1.0mM (0.2924g/L)

Store at 4°C

## Citrate Synthase

Reagent mix per mL

Assay buffer 500µL

Acetyl - CoA 20

Ellman's reagent 20

Triton X-100 10% 10

H<sub>2</sub>O 500

Heat reagent mix to 30°C in water bath.

per well

Homogenate variable (~5-10µL)

Oxaloacetic acid 7µL (0 for blank)

Reagent Mix 200µL

Measure increase in absorbance at 412 nm at 30°C

## Reagents

-Assay buffer- Tris 200mM (24.24g/L) pH 7.4 with HCl. Store at 4°C

-Acetyl CoA (Roche 101907)- 8.3mg/1.0mL (5.5mL per 50mg bottle). Prepare in Milli Q H<sub>2</sub>O. Store at -20°C indefinitely

-Ellman's (DTNB) (Sigma D8130)- 4.0mg /1.0mL. Prepare 40mg in 10mL 95% ETOH. Store at -20°C indefinitely

-Oxaloacetic Acid (Sigma O-4126) 1.3mg/1.0mL. Prepare 10mL and pH to 7.4 with KOH. Aliquot and store at -20°C for one month.

-Triton X-100 10% in H<sub>2</sub>O v/v. Store at 4°C or -20°C indefinitely.

## Cytochrome C Oxidase

Reagent Mix per mL

Assay buffer	500 $\mu$ L
Cytochrome C	variable
Dodecylmaltoside 10%	10 $\mu$ L
H <sub>2</sub> O	500 $\mu$ L

Heat reagent mix to 30°C in water bath.

per well

Homogenate	variable (~5-10 $\mu$ L)
Reagent Mix	200 $\mu$ L

Measure decrease in absorbance at 550nm at 30°C

Reagents

-Assay buffer	Potassium Phosphate 100mM pH7.0
	K <sub>2</sub> HPO <sub>4</sub> 10.71g/L H <sub>2</sub> O
	KH <sub>2</sub> PO <sub>4</sub> 5.24g/L H <sub>2</sub> O
	Store at 4°C

Reduced Cytochrome C- cytochrome c (SigmaC2506) 123.5mg/mL in 50mM potassium phosphate buffer pH 7.0 . Add 7.5 mg/mL sodium dithionite and incubate in the dark on ice for 5 minutes. Purify on a 10K MWC0 spin filter (Millipore UFV2BGC10) by spinning at 2000g for ~one hour, adding more buffer every 15 minutes and mixing with the cytochrome c. Check the reduction of the cytochrome c by measuring the ratio of absorbance at 550/565nm. It should be >6. To calculate the concentration of reduced cytochrome c (19.6)

n-Dodecylmaltoside- 10% w/v in H<sub>2</sub>O. Store at 4°C or -20°C indefinitely.