

Northern Blotting

Solutions

RNase-Free water

Prepare DEPC 0.1%(v/v) in Milli-Q water. Mix well, autoclave.

10X MOPS/EDTA

0.5M MOPS pH 7.0

0.01M EDTA pH 7.5 – 8.0

For 500mL: 53.325g MOPS
10mL 0.5M EDTA pH 8.0
Bring up to 500mL with Milli-Q
pH to 7.0 with NaOH
Filter through 0.2 μ M filter
Store wrapped in aluminum foil

Formaldehyde/Formamide

89 μ L 37% formaldehyde
250 μ L deionized formamide

Buffer A

294 μ L 10X MOPS/EDTA
106 μ L RNase-free water

Gel Loading Buffer

322 μ L Buffer A
5mg xylene cyanol FF
5mg bromocresol green
400mg sucrose

Mix well, add:

178 μ L 37% formaldehyde
500 μ L formamide

Sample Preparation

1. Speed vac (on "medium" drying rate) 10µg total RNA, or 1-2µg poly A RNA. The lids of the tubes should be open, and their mouths covered with perforated parafilm.
2. Add 5µL Buffer A and resuspend RNA.
3. Add 12µL formaldehyde/formamide solution and mix.
4. Incubate samples at 68°C for 10 minutes and place on ice.
5. Add 4µL gel loading buffer and mix.

Gel Preparation and Electrophoresis

1. Wash tray, comb and electrophoresis apparatus with 10% SDS and rinse thoroughly with Milli-Q. Alternatively, soak them in RNase free water overnight, or *at least* while performing the subsequent steps.
2. Prepare a 1% agarose gel in 1X MOPS/EDTA buffer and cool to 60°C.
3. In a fume hood add 9mL 37% formaldehyde to 41 ml of the cooled agarose. Mix and cast gel. Allow gel to polymerize for at least 45 minutes.
4. Pre-run gel at 40 volts for 45 minutes in 1X MOPS/EDTA buffer.
5. Load samples and run gel at 40 volts for small gels (80 volts for large gels).
6. Allow gel to run until bromophenol blue dye reaches ~ 1cm from the bottom.

Transfer

1. After electrophoresis, wash gel in 20X SSC for 20 minutes.
2. Pre-wet membrane in Milli Q water or in 2X SSC.
3. Transfer gel using Vacugene apparatus. The order of components, starting from the bottom, is: screen, membrane, mask, gel. The membrane should be ~ 1 cm larger than the opening of the mask all around. Fill with 20X SSC (10X SSC works equally well) until gel is covered. Transfer at 70—80 cm/H₂O pressure for 2 hours.
4. After transfer, rinse membrane in 2X SSC and crosslink ("Auto Cross Link" function of the Stratalinker).

Prehybridization

Pre-hyb solution (final concentrations):

Volume of stock solutions required:

55% formamide

5.5 ml formamide

5X SSC

2.5 ml 20X SSC

1x Denhardtts

200 µl 50X Denhardtts

50mM phosphate buffer, pH 7.0

2.5mL 0.2M phosphate buffer, pH7.0

250µg/mL ssDNA (fish sperm DNA)

250 µl 10 mg/ml ss DNA

-Alternatively, 500µL 50X Denhardtts and 1.5mL 0.2M phosphate buffer, pH 7.0 (other components the same) will work equally well.

-For 100mL 0.2M phosphate buffer (pH 7.0), mix 57.7mL 1M Na₂HPO₄ and 42.3mL 1M NaH₂PO₄

-Incubate the membrane in the pre-hyb solution for *at least* 4 hours at 42°C (overnight, if preferred).

Hybridization

Hyb solution (final concentrations):

Volume of stock solutions required:

55% formamide	3.3mL formamide
5X SSC	1.5mL 20X SSC
1x Denhardt's	120µL 50X Denhardt's
20mM phosphate buffer, pH 7.0	600µL 0.2 M phosphate buffer, pH 7.0
10% dextran sulphate	1.2mL 50% dextran sulphate
250µg/mL ssDNA (salmon sperm DNA)	150µL 10mg/mL ssDNA

-For labeling and counting of the probe, see protocol for Southern blotting.

-Prepare hyb mixture without the probe. Place the membrane in the container; add the hyb solution, and finally, the probe. Incubate overnight at 42°C.

Washing

Two alternative methods (work equally well):

Method 1

-Remove membrane from hyb solution and rinse 3x at room temperature with 2X SSC+0.1% SDS solution.

-Wash the membrane 2x 30 minutes, in 2X SSC+0.1% SDS, at 42°C while shaking.

Method 2

-Remove membrane from hyb solution and rinse 3x at room temperature with 2X SSC+0.1% SDS solution.

-Wash the membrane 4x 5 minutes in 2X SSC+0.1% SDS at room temperature while shaking.

-Wash the membrane 2x 15 minutes in 0.1X SSC+0.1% SDS at 50°C while shaking.

-Regardless of the method used, check counts with Geiger in between the longer washes, especially if the probe is not particularly "hot". If the membrane is still "hot" overall, continue washing.

Otherwise, stop and expose the membrane to a phosphorimager screen for at least 4 hours.

Stripping

-Place membrane in 0.1X SSC+0.1% SDS, 2x 30 minutes, at 85°C with vigorous shaking.